

TELOS Arthroscopy Leg Holder


The TELOS Leg Holder is a low profile, easy to use holder providing simple velcro release and fixation of the leg during arthroscopic procedures. Telos Leg Holder allows use of standard tourniquet cuffs and can be used with any standard O.R. Table system. It easily adjusts to fit any size leg with the simple velcro securing strap and is radiolucent in design.



RLH630


OTHER TELOS SUPPORTS:

Anesthesia Arm Board




RAS610

Arm Support



RAB620

Radiolucent
-Hand, -Wrist Board



RHB640

Quality Orthopedic Products

TELOS Medical / Austin & Associates
212 Copperwood Court, Millersville, MD 21108
PH (410) 903-9038 • FX (410) 877-0544
Email: telousa@aol.com | www.telosmedical.com